

Screening Drugs and Toxic Compounds with LC-MS/MS: An Alternative to LC-UV for Research Toxicology Labs

Jordan Velardo¹, Monique Manchon¹, Bénédicte Duret², Dennis Nagtalon³, Marta Kozak³;

¹Laboratory of Toxicology, Lyon Sud Hospital, Pierre-Bénite, France; ²Thermo Fisher Scientific, Les Ulis, France; ³Thermo Fisher Scientific, San Jose, CA, USA

Key Words

- ToxSpec Analyzer
- ToxID software
- LXQ Linear Ion Trap
- Accela UHPLC System

Introduction

Screening for drugs of abuse and other toxic compounds in biological samples has quickly become a routine assay conducted in many research toxicology laboratories. The main challenge is to get rapid and accurate results amidst the generally large number of potential analytes to be identified within complex biological matrices. One of the techniques widely used in this area is high pressure liquid chromatography (HPLC) combined with photo diode array detection (DAD) or ultra-violet (UV) detection. The most popular LC-UV platform has been the Bio-Rad® REMEDI™ HS drug profiling system. When this platform was recently discontinued, a significant technological gap became apparent. Now this gap is rapidly being filled by newer, more effective high pressure liquid chromatography - mass spectrometry (HPLC-MS) technologies.

Here we present the workflow and results obtained by using the Thermo Scientific ToxSpec Analyzer, a new UHPLC-MS system based on ultra-high pressure liquid chromatography and linear ion trap mass spectrometry technology.

Goal

Evaluate an LC-MS/MS method for screening and semi-quantitation of drugs and toxic compounds in serum and urine matrices to determine if this approach can provide an alternative to REMEDI technology for research toxicology.

Experimental

The ToxSpec™ Analyzer combines hardware, software, and screening methodologies designed to significantly simplify and improve the screening assay workflow. LC-MS² data is acquired by using a pre-

configured instrument method, and the data is automatically processed, post-acquisition, by Thermo Scientific ToxID automated drug screening software.

The LC-MS screening was performed on Thermo Scientific instrumentation including an LXQ™ linear ion trap mass spectrometer coupled to an Accela™ UHPLC system using a polarity-switching and scan-dependent MS/MS experiment (Figure 1). The MS² spectra generated were processed through ToxID™ software. Using a novel screening algorithm, the software program identifies target analytes through a MS² library search against a large spectral library of known analytes as well as expected retention times. Semi-quantitative data results can also be generated concurrently from the MS² spectral intensity ratios between the target analyte and the corresponding internal standard.




Figure 1: Polarity-switching and scan-dependent MS/MS experiment

The ToxSpec Analyzer includes a diverse and easily-expandable MS/MS library of 300 compounds that it screens using a single pre-configured method. In our laboratory, we have expanded the library by more than 50 entries to date.

Sample preparation

The extraction procedure was performed by using liquid/liquid extraction (LLE) with Toxi-Tube A® (Varian, les Ulis, France). Details of the procedure are described below.

- Vortex the Toxi-Tube A for 10 seconds.
- Add 1 mL of serum or urine into the Toxi-Tube A.
- Add 200 µL of a solution of internal standard [haloperidol-d4, chlorpromazine-d3, and prazepam-d5 at the following concentrations: 100 ng/mL, 1 µg/mL and 100 ng/mL, respectively, in 70/30 of A/B (A: water containing 10 mM ammonium acetate and 0.1% formic acid; B: acetonitrile containing 0.1% formic acid)].
- Add 5 mL of water.
- Vortex for 10 seconds.
- Mix for 5 minutes.
- Centrifuge for 5 minutes at 2700 rpm.
- Transfer the upper layer to a tube and evaporate to dryness at 40 °C.
- Reconstitute the sample in 200 µL of 70/30 of A/B.

HPLC Conditions

Chromatographic analyses were performed using the Thermo Scientific Accela UHPLC system. The chromatographic conditions were as follows:

Column:	Thermo Scientific Hypersil GOLD PFP 5 µm, 150 x 2.1 mm		
Flow rate:	0.2 mL/min		
Mobile phase:	A: water containing 10 mM ammonium acetate and 0.1% formic acid; B: acetonitrile containing 0.1% formic acid		
Injection volume:	10 µL		
Gradient:	T (min)	A (%)	B (%)
	0.0	95	5
	5.0	55	45
	18.0	30	70
	20.0	5	95
	27.0	5	95
	27.1	95	5
	32.0	95	5

Mass Spectrometry Conditions

MS analysis was carried out on our LXQ linear ion trap mass spectrometer with an electrospray ionization (ESI) source. The MS conditions were as follows:

Ion polarity:	Polarity-switching scan-dependent experiment
Spray voltage:	5000 V
Sheath gas (N ₂) pressure:	30 (arbitrary units)
Auxiliary gas (N ₂) pressure:	8 (arbitrary units)
Capillary temperature:	275 °C
Microscan:	1
Wideband Activation™:	Activated
Stepped normalized collision energy:	35% ± 10%

Results and Discussion

More than 150 real laboratory samples (serum and urine) have been analyzed. Table 1 reports some of the data obtained from both the REMEDi HS LC/UV system and the ToxSpec Analyzer UHPLC/MS system. Among the 12 samples reported here, 22 compounds have been identified using both the REMEDi HS and the ToxSpec Analyzer. Notably however, the ToxSpec Analyzer system identified 24 additional compounds that were not detected with the REMEDi HS due in most cases to a lack of sensitivity, specificity, and coelution capability.

The ToxSpec Analyzer also provided a better response for some classes of compounds, like benzodiazepines. With the REMEDi HS system, the retention time for this class of compounds was close to the dead volume of the column. For that reason, the signals that interfered with matrix components were rather difficult to identify. It was also observed that haloperidol (sample #5) and paroxetine (sample #10) gave a much better signal with the ToxSpec Analyzer.

Sample #	Compounds identified using ToxSpec Analyzer	Compounds identified using REMEDi HS
1	Acetaminophen Nortriptyline Amitriptyline Oxazepam	Not detected Not detected Amitriptyline Not detected
2	Nordiazepam Alprazolam Cyamemazine	Nordiazepam Not detected Cyamemazine
3	Acetaminophen Nordiazepam Venlafaxine Oxazepam Alprazolam	Not detected Nordiazepam Venlafaxine Oxazepam Not detected
4	Nordiazepam Diazepam Oxazepam Temazepam Levomepromazine Zopiclone	Not detected Diazepam Not detected Not detected Levomepromazine Zopiclone
5	Oxazepam Clomipramine Quinidine Haloperidol Clonazepam	Not detected Clomipramine Quinine Not detected Not detected
6	Acetaminophen Bisoprolol	Not detected Bisoprolol
7	Venlafaxine Risperidone	Venlafaxine Not detected
8	Quinine Hydromorphone Morphine	Quinine Hydromorphone Morphine
9	Lidocaine Nortriptyline Mirtazapine Amitriptyline Cyamemazine Levomepromazine Zopiclone	Not detected Not detected Not detected Amitriptyline Cyamemazine Levomepromazine Not detected
10	Bromazepam Paroxetine	Bromazepam Not detected
11	Sertraline Hydrocortisone	Not detected Not detected
12	Acetaminophen Alprazolam Prednisolone Hydroxyzine Fexofenadine	Not detected Alprazolam Not detected Hydroxyzine Not detected
TOTAL	46 Molecules	22 molecules

Table 1: List of psychoactive molecules identified in real laboratory samples using the ToxSpec Analyzer compared to the REMEDi HS system

Our aim was to quickly and confidently identify toxic compounds in the samples by spectral library searching while performing a semi-quantification calculation for identified compounds. To perform the semi-quantification, a response factor that correlated the intensity of the MS² spectra to a concentration was calculated for each molecule present in the library using internal standards. The semi-quantification result was automatically calculated using ToxID software. An example of the automatically-generated report can be seen in Figure 2. The report includes a list of compounds identified in a real laboratory sample and their respective calculated concentrations.


One important aspect of this method is the ability to reprocess data retrospectively from the MS spectra. The ToxID report is based on MS² spectra library searching. This means that if the entry corresponding to the compound is not currently available in the library, ToxID will not be able to identify the analyte. However, as data are acquired in MS mode, it is possible to reprocess the MS trace and check that all major ions have been identified by ToxID. If not, it is then possible to re-inject the sample and perform MS² acquisition on specific ions.

Conclusion

The ToxSpec Analyzer is a good replacement for the REMEDi HS system in research toxicology laboratories because it offers increased sensitivity, greater specificity, and lower cost-per sample analysis.

Centre Hospitalier Lyon Sud Summary Report

Raw File Name: C:\Xcalibur\Validation\REMED\Real samples\090107\375489.RAW
 Config File Name: C:\Xcalibur\Validation\REMED\CSV fi...LXQ_GUS_config_30min_TOXID Semi-quant.csv
 Sample Name: Laboratory: Laboratoire de Toxicologie
 Acquisition Start Time: 1/7/2009 1:23:18 PM
 Screening Conditions: Based on Full MS2 scans. m/z window (amu): 0.50. RT window (min): 2.00. MS2 Search libraries: Tox_Library_LXQ. Use full MS scan to confirm.



Peak Number	Compound Name	Code	SI	RSI	m/z	Expected RT	Actual RT	Concentration mg/l	Library Name
1	Nordiazepam	p	700	708	271.10	8.60	9.38	1.05	Tox_Library_LXQ
2	Alprazolam	p	811	814	309.10	9.10	9.61	0.84	Tox_Library_LXQ
3	Cyamemazine	p	818	819	324.20	15.00	14.79	0.25	Tox_Library_LXQ
4	Prazepam-D5	i	890	894	330.20	11.70	12.25	0.87	Tox_Library_LXQ
5	Hydrocortisone	p	863	865	363.20	7.50	7.72	0.11	Tox_Library_LXQ
6	Haloperidol-D4	i	845	861	380.40	14.90	14.60	1.00	Tox_Library_LXQ




Figure 2: ToxID report – short summary style

Legal Notices

©2009 Thermo Fisher Scientific Inc. All rights reserved. REMEDi is a trademark of Bio-Rad Laboratories. Toxi-Tube A is a registered trademark of Varian, Inc. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. LC-MS Systems are for Research Use Only. Not to be used in Diagnostic Procedures.

View additional Thermo Scientific LC/MS application notes at: www.thermo.com/appnotes

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

- Africa-Other**
+27 11 570 1840
- Australia**
+61 2 8844 9500
- Austria**
+43 1 333 50 34 0
- Belgium**
+32 2 482 30 30
- Canada**
+1 800 530 8447
- China**
+86 10 8419 3588
- Denmark**
+45 70 23 62 60
- Europe-Other**
+43 1 333 50 34 0
- Finland/Norway/Sweden**
+46 8 556 468 00
- France**
+33 1 60 92 48 00
- Germany**
+49 6103 408 1014
- India**
+91 22 6742 9434
- Italy**
+39 02 950 591
- Japan**
+81 45 453 9100
- Latin America**
+1 608 276 5659
- Middle East**
+43 1 333 50 34 0
- Netherlands**
+31 76 579 55 55
- South Africa**
+27 11 570 1840
- Spain**
+34 914 845 965
- Switzerland**
+41 61 716 77 00
- UK**
+44 1442 233555
- USA**
+1 800 532 4752

www.thermo.com



Thermo Fisher Scientific, San Jose, CA USA is ISO Certified.

AN63034_E 07/09S