



Buffers for Life Science Research

CellPURE™ Biological Buffers | JustPURE™ "Good" Buffers | Core Essential Buffers

Quality.
Consistency.
Selection.



- Ultrapure zwitterionic biological buffers
- Optimized for cell/tissue culture work
- Wide applicability to molecular biology and biochemical studies
- Manufactured under strict quality control guidelines to ensure performance and reliability

Fisher BioReagents buffers meet the needs of every budget and scale.

- Economical powders come in various package sizes
- Concentrated stock solutions provide convenience in one step of easy dilution
- Ready-to-use solutions offer the biggest time savings of all

Obtaining optimal results in your research requires careful selection of reagents. When your experiments require exact buffering conditions, you can depend on the reliability of Fisher BioReagents buffers. All Fisher BioReagents buffers are manufactured from high-quality raw materials under ISO 9001:2008-certified manufacturing and testing processes.

The purpose of a buffer in a biological system is to maintain intracellular and extracellular pH within the physiological range and resist changes in pH due to the presence of internal and external influences. CellPURE Biological Buffers from Fisher BioReagents are ultrapure zwitterionic buffers that possess both positive and negative charges. First described by Good and co-workers in 1966, these organic buffers have pKa values at or near physiological pH and minimal interference with biological processes. CellPURE Biological Buffers are ideal for cell cultivation, isolation of cells, enzyme assays, and other biochemical applications.

(Reference: Good, N.E., et al. (1966) Hydrogen Ion Buffers for Biological Research. *Biochemistry* 5(2): 467-477)

Advantages

- Cell culture tested
- Analyzed for the absence of nucleases and proteases
- Tested for endotoxin and bioburden levels
- pK_a values mostly independent of temperature and concentration
- High water solubility
- Minimal permeability to biological membranes

Applications

- Tissue culture media and maintenance of cell lines
- Enzyme assays
- Electrophoresis of RNA
- Transfection of mammalian cells

Properties of CellPURE Biological Buffers

Cat. No.	Product Description	Formula	Molecular Weight	pK _a @ 25° C	Useful pH Range	Cell Culture Tested	Endotoxin Assay	Nuclease and Protease Tested	Size
BP2941-100	BIS-TRIS	C ₈ H ₁₉ NO ₅	209.24	6.5	5.8 to 7.2	X	X	X	100g
BP2943-100	BIS-TRIS propane	C ₁₁ H ₂₆ N ₂ O ₆	282.33	6.8†	6.3 to 9.5	X	X	X	100g
BP2947-100	BES	C ₆ H ₁₅ NO ₅ S	213.25	7.1	6.4 to 7.8	X	X	X	100g
BP2936-100	MOPS	C ₇ H ₁₅ NO ₄ S	209.26	7.2	6.5 to 7.9	X	X	X	100g
BP2946-25	MOPS sodium salt	C ₇ H ₁₄ NNaO ₄ S	231.25	7.2	6.5 to 7.9	X	X	X	25g
BP2937-100	HEPES	C ₈ H ₁₈ N ₂ O ₄ S	238.30	7.5	6.8 to 8.2	X	X	X	100g
BP2939-100	HEPES Sodium Salt	C ₈ H ₁₇ N ₂ NaO ₄ S	260.29	7.5	6.8 to 8.2	X	X	X	100g
BP2945-100	TES	C ₆ H ₁₅ NO ₆ S	229.25	7.5	6.8 to 8.2	X	X	X	100g

† pK_a = 9.0 for the second dissociation stage

Technical Tip

The Water Makes a Difference. CellPURE Biological Buffers are of the highest quality (purity ≥ 99%). They are tested for low heavy metal content and the absence of endotoxins, nucleases, and proteases. It is important to use only water of the highest quality to prepare the buffer solutions. Water that stands too long in pipes or water produced by a system needing routine maintenance increases the risk for contamination of buffer solutions. Fisher BioReagents provides several grades of high purity water suitable for the preparation of high quality buffer solutions.



Zwitterionic buffers possessing both positive and negative charges were described in 1966 by Good and co-workers as being suitable for work with biological molecules. Popularly known today as "Good" buffers, these organic buffers have several advantages compared to traditional inorganic buffering systems such as phosphate, borate, and bicarbonate. "Good" buffers approach the "ideal" buffer state by having pKa values at or near physiological pH and minimal interference with biological processes. JustPURE "Good" buffers from Fisher BioReagents have very high purity (assay > 99%) and only trace amounts of metal ions, which makes them useful for applications requiring tight control of elemental content.

(Reference: Good, N.E., et al. (1966) Hydrogen Ion Buffers for Biological Research. *Biochemistry* 5(2): 467-477)

- Ultrapure zwitterionic buffers
- Optimized for research in cellular and molecular biology
- Low interference with biological reactions

- pKa values mostly independent of temperature and concentration
- High water solubility
- Chemically and enzymatically stable
- Minimal permeability to biological membranes
- Limited effect on biochemical reactions
- Minimal absorption in spectral range 240 to 700nm


Applications


- High efficiency transfection of mammalian cells
- Gel electrophoresis of RNA
- Protein isolation applications
- Cell cultures and enzyme assays
- Bioanalytical methods such as IEF, 2-D electrophoresis, SDS-PAGE

Advantages Properties of JustPURE Buffers

Cat. No.	Product Description	CAS No.	Formula	Molecular Weight	pK _a @ 25° C	Useful pH Range	Form	Size
BP2920-250	MES hydrate	4432-31-9	C ₆ H ₁₃ NO ₄ S·xH ₂ O	195.24 (anhydrous)	6.1	5.5 to 6.7	White powder	250g
BP2935-100	ACES	7365-82-4	C ₄ H ₁₀ N ₂ O ₄ S	182.20	6.8	6.1 to 7.5	White powder	100g
BP2924-50	PIPES	5625-37-6	C ₈ H ₁₈ N ₂ O ₆ S ₂	302.37	6.8	6.1 to 7.5	White powder	50g
BP2929-25	BIS-TRIS propane	64431-96-5	C ₁₁ H ₂₆ N ₂ O ₆	282.33	6.8†	6.3 to 9.5	White powder	25g
BP2925-100	MOPS	1132-61-2	C ₇ H ₁₅ NO ₄ S	209.26	7.2	6.5 to 7.9	White powder	100g
BP2921-50	HEPES	7365-45-9	C ₈ H ₁₈ N ₂ O ₄ S	238.30	7.5	6.8 to 8.2	White powder	50g
BP2933-100	EPPS	16052-06-5	C ₉ H ₂₀ N ₂ O ₄ S	252.33	8.0	7.3 to 8.7	White powder	100g
BP2922-100	Tricine	5704-04-1	C ₆ H ₁₃ NO ₅	179.17	8.1	7.4 to 8.8	White powder	100g
BP2930-50	Gly-Gly	556-50-3	C ₄ H ₈ N ₂ O ₃	132.12	8.2	7.5 to 8.9	White powder	50g
BP2923-100	TAPS	29915-38-6	C ₇ H ₁₇ NO ₆ S	243.28	8.4	7.7 to 9.1	White powder	100g
BP2934-100	CHES	103-47-9	C ₈ H ₁₇ NO ₃ S	207.29	9.3	8.6 to 10.0	White powder	100g
BP2928-100	CAPS	1135-40-6	C ₉ H ₁₉ NO ₃ S	221.32	10.4	9.7 to 11.1	White powder	100g

† pK_a = 9.0 for the second dissociation stage





Choose Fisher BioReagents Buffers for:

- QUALITY:** tight specifications
- CONSISTENCY:** lot-to-lot uniformity
- SELECTION:** powders, concentrated stock solutions, or ready-to-use liquids
- PACKAGING:** state-of-the-art containers designed for safety and utility
- CONVENIENCE:** pre-qualified for a variety of applications
- ECONOMY:** configurations to meet all budgets
- SCALE:** from bench to batch sizes

Fisher BioReagents buffers meet the needs of every budget and scale.

- Economical powders come in various package sizes
- Concentrated stock solutions provide convenience in one step of easy dilution
- Ready-to-use solutions offer the biggest time savings of all

Suitable for Life Science applications:

- Electrophoresis of nucleic acids and proteins
- Isolation and purification of biomolecules
- Cell and tissue culture
- Guaranteed DNase-, RNase-, and protease-free

Obtaining optimal results in your research requires careful selection of reagents. When your experiments require exact buffering conditions, you can depend on the reliability of Fisher BioReagents buffers. All Fisher BioReagents buffers are manufactured from high-quality raw materials under ISO 9001-2008-certified manufacturing and testing processes. In addition, buffer sterilization is done by 0.2µm filtration and/or autoclaving.

Description	Size	Cat. No.
Buffers for Protein Electrophoresis		
TG Tris-Glycine		
10X	1L	BP1306-1
10X	4L	BP1306-4
10X	1L**	BP1307-1
TGS Tris-Glycine-SDS		
10X	1L	BP1341-1
10X	4L	BP1341-4
5X	1L*	BP1398-92
10X	1L**	BP1342-1
SDS Sodium Dodecyl Sulfate		
10%	200mL	BP2436-200
10%	1L	BP2436-1
20%	200mL	BP1311-200
20%	1L	BP1311-1
White Powder	100g	BP166-100
White Powder	500g	BP166-500
White Powder	5kg	BP166-5
Buffers for Nucleic Acid Electrophoresis		
TBE Tris-Borate-EDTA		
1X	1L	BP2430-1
1X	4L	BP2430-4
1X	20L	BP2430-20
5X	1L*	BP1396-86
10X	1L	BP1333-1
10X	4L	BP1333-4
10X	20L	BP1333-20
10X	1L**	BP1334-1
TAE Tris-Acetate-EDTA		
1X	4L	BP2434-4
1X	20L	BP2434-20
10X	500mL	BP1335-500
10X	1L	BP1335-1
10X	4L	BP1335-4
10X	20L	BP1335-20
25X	1L	BP1330-1
50X	500mL	BP1332-500
50X	1L	BP1332-1
50X	4L	BP1332-4
50X	20L	BP1332-20
25X	1L**	BP1331-1
MOPS Buffer		
Powder	100g	BP308-100
Powder	500g	BP308-500
10X Buffer Solution	500mL	BP2900-500
10X Buffer Solution	1L	BP2900-1
Buffers for Hybridization of Nucleic Acids		
SSPE Saline-Sodium Phosphate-EDTA		
20X	1L	BP1328-1
20X	4L	BP1328-4
20X	20L	BP1328-20
SSC Saline-Sodium Citrate		
20X	1L	BP1325-1
20X	4L	BP1325-4
20X	20L	BP1325-20

Choose Fisher BioReagents® buffers for:

QUALITY: tight specifications

CONSISTENCY: lot-to-lot uniformity

SELECTION: powders, concentrated stock solutions, or ready-to-use liquids

PACKAGING: state-of-the-art containers designed for safety and utility

CONVENIENCE: pre-qualified for a variety of applications

ECONOMY: size configurations to meet all budgets

SCALE: from bench to batch sizes

Description	Size	Cat. No.
Buffers for Life Science Applications		
EDTA Ethylenediamine Tetraacetic Acid		
0.5M (pH 8.0)	100mL	BP2482-100
0.5M (pH 8.0)	500mL	BP2482-500
0.5M (pH 8.0)	1L	BP2482-1
0.5M (pH 8.0)	20L	BP2482-20
0.5M (pH 8.0), DEPC	100mL	BP2483-100
0.5M (pH 8.0), DEPC	500mL	BP2483-500
0.5M (pH 8.0), DEPC	1L	BP2483-1
Powder	500g	BP118-500
HEPES		
1.0M (pH 7.3)	100mL	BP299-100
1.0M (pH 7.3)	500mL	BP299-500
1.0M (pH 7.3)	1L	BP299-1
Crystals	100g	BP310-100
Crystals	500g	BP310-500
Crystals	1kg	BP310-1
Crystals	5kg	BP310-5
PBS Phosphate Buffered Saline		
Tablets	100 tablets	BP2944-100
1X	4L	BP2438-4
1X	20L	BP2438-20
10X	500mL	BP399-500
10X	1L	BP399-1
CellIPURE™ PBS 10X	4L	BP2940-4
10X	4L	BP399-4
10X	20L	BP399-20
1X Powder Concentrate	10L	BP661-10
1X Powder Concentrate	50L	BP661-50
Powder with Tween® 20	10 pouches	BP2938-10
Powder with BSA	10 pouches	BP2942-10
10X Powder Concentrate	2x1L	BP665-1
STE Sodium Chloride-Tris-EDTA		
1X (pH 8.0)	1L	BP2478-1
10X (pH 8.0)	1L	BP2479-1
TBS Tris-Buffered Saline		
1X (pH 7.4)	1L	BP2472-1
10X (pH 7.4)	100mL	BP2471-100
10X (pH 7.4)	1L	BP2471-1

Description	Size	Cat. No.
TE Tris-EDTA		
1X (pH 7.4)	100mL	BP2476-100
1X (pH 7.4)	500mL	BP2476-500
1X (pH 7.4)	1L	BP2476-1
1X (pH 8.0)	100mL	BP2473-100
1X (pH 8.0)	500mL	BP2473-500
1X (pH 8.0)	1L	BP2473-1
1X (pH 7.6)	100mL	BP2474-100
1X (pH 7.6)	500mL	BP2474-500
1X (pH 7.6)	1L	BP2474-1
10X (pH 7.4)	1L	BP2477-1
10X (pH 7.6)	100mL	BP2475-100
10X (pH 7.6)	500mL	BP2475-500
10X (pH 7.6)	1L	BP2475-1
100X (pH 8.0)	1L	BP1338-1
100X (pH 8.0)	4L	BP1338-4
100X	1L**	BP1339-1
Tris Buffer		
0.3M	500mL	BP1761-500
0.3M	1L	BP1761-1
2.0M	100mL	BP1759-100
2.0M	500mL	BP1759-500
Tris Hydrochloride		
1.0M (pH 7.0)	100mL	BP1756-100
1.0M (pH 7.0)	500mL	BP1756-500
1.0M (pH 7.5)	100mL	BP1757-100
1.0M (pH 7.5)	500mL	BP1757-500
1.0M (pH 8.0)	100mL	BP1758-100
1.0M (pH 8.0)	500mL	BP1758-500
Solid	500g	BP153-500
Solid	1kg	BP153-1
Tris Base		
Crystals	500g	BP152-500
Crystals	1kg	BP152-1
Crystals	5kg	BP152-5
Crystals	10kg	BP152-10
Crystals	25kg	BP152-25
Water		
Biotech Grade	4L	BP2485-4
	20L	BP2485-20
Nuclease-Free, Sterile	50mL	BP2484-50
	100mL	BP2484-100
DNA Grade, Sterile	1L	BP2470-1
RNA Grade, DEPC-treated, Sterile	1L	BP561-1

Quality. Consistency. Selection.

* Pre-weighed powder in poly bottle. Dissolve in water.

** Pre-weighed powder in foil pack. Dissolve in water.

Related Core BioReagents

Ethanol,

Molecular Biology Grade

Cat. No.	Size
BP2818-100	100mL
BP2818-500	500mL
BP2818-4	4L



Isopropanol,

Molecular Biology Grade

Cat. No.	Size
BP2618-500	500mL
BP2618-1	1L
BP2618-212	2.5L
BP2618-4 4L	

Water (0.1µm filtered),

Molecular Biology Grade

Cat. No.	Size
BP2819-100	100mL
BP2819-1	1L
BP2819-4	4L
BP2819-10	10L
BP2819-20	20L





To place an order, contact your local Fisher Scientific office.



Fisher BioReagents™

©2014 Thermo Fisher Scientific Inc. All rights reserved.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

AMERICAS

Canada
Fisher Scientific Canada
www.fishersci.ca
Tel: 800-234-7437
Fax: 800-463-2996

Latin America
Fisher Scientific Global Export,
Latin America
www.fishersci.com
Tel: 770-871-4725
Fax: 770-871-4726

United States
Fisher HealthCare
www.fisherhealthcare.com
Tel: 800-640-0640
Fax: 800-290-0290

Fisher Scientific
www.fishersci.com
Tel: 800-766-7000
Fax: 800-926-1166

EUROPE

Austria
info.austria@thermofisher.com
www.at.fishersci.com
Tel: 0800 20 88 40
Fax: 0800 20 66 90

Belgium
be.fisher@thermofisher.com
www.be.fishersci.com
Tel: 056 260 260
Fax: 056 260 270

Czech Republic
info.cz@thermofisher.com
www.thermofisher.cz
Tel: 466 798 230
Fax: 466 435 008

Denmark
kundeservice@thermofisher.com
www.dk.fishersci.com
Tel: 70 27 99 20
Fax: 70 27 99 29

Finland
fisher.fi@thermofisher.com
www.fi.fishersci.com
Tel: 09-802 76 280
Fax: 09-802 76 235

France
fr.fisher@thermofisher.com
www.fr.fishersci.com
Tel: 03 88 67 14 14
Fax: 03 88 67 11 68

Germany
info.germany@thermofisher.com
www.de.fishersci.com
Tel: 2304 932-5
Fax: 2304 932-950

Ireland
fsie.sales@thermofisher.com
www.ie.fishersci.com
Tel: 01 885 5854
Fax: 01 899 1855

Italy
it.fisher@thermofisher.com
www.it.fishersci.com
Tel: 02 950 59 478
Fax: 02 950 59 479

The Netherlands
nl.info@thermofisher.com
www.nl.fishersci.com
Tel: 020 487 70 00
Fax: 020 487 70 70

Norway
fisher.no@thermofisher.com
www.no.fishersci.com
Tel: 22 95 59 59
Fax: 22 95 59 40

Portugal
pt.fisher@thermofisher.com
www.pt.fishersci.com
Tel: 21 425 33 50
Fax: 21 425 33 51

Spain
es.fisher@thermofisher.com
www.es.fishersci.com
Tel: 902 239 303
Fax: 902 239 404

Sweden
fisher.se@thermofisher.com
www.se.fishersci.com
Tel: 31 352 32 00
Fax: 31 352 32 50

Switzerland
info.ch@thermofisher.com
www.ch.fishersci.com
Tel: 056 618 41 11
Fax: 056 618 41 41

United Kingdom
fsuk.sales@thermofisher.com
www.fisher.co.uk
Tel: 01509 555500
Fax: 01509 555111

ASIA

India
qfc.customercare@thermofisher.com
www.fishersci.in
Tel: +91-22-6680 3000
Fax: +91-22-6680 3001

Korea
sales.fsk@thermofisher.com
www.fishersci.co.kr
Tel: 02-3420-8700
Fax: 02-3420-8710

Malaysia
enquiry.my@thermofisher.com
www.fishersci.com.my
Tel: (603) 51228888
Fax: (603) 5121 8899

Singapore
enquiry.sg@thermofisher.com
www.fishersci.com.sg
Tel: (65) 6873 6006
Fax: (65) 6873 5005

OCEANIA

Australia
AUinfo@thermofisher.com
www.thermofisher.com.au
Tel: 1300-735-292
Fax: 1800-067-639

New Zealand
NZinfo@thermofisher.com
www.thermofisher.co.nz
Tel: 0800-933-966
Fax: +64 9 980 6788

