

Designed for performance Optimised for comfort

Deliver high quality routine results comfortably and efficiently with the Thermo Scientific™ HM525 NX Cryostat, a high performance routine cryostat with intuitive software and touchscreen for simple, efficient operation. The HM525 NX Cryostat combines the proven operation and reliability of the HM525 cryostat and the new ergonomics and intuitive interface of the Thermo Scientific NX series cryostats.

Exceptional sectioning quality and reliability

Reliable stepper motor technology delivers reproducible section thickness

User friendly touch screen interface

Icon driven interface allows intuitive operation

Enhanced ergonomics

Light-touch handwheel requires minimal force to operate

Operator safety

Mechanical handwheel lock completely immobilises the specimen head

27 cryobar storage positions including 4 quick freeze positions

Allows for greater workflow flexibility and improved specimen management

Optional UV disinfection



Ordering Information

Code	Product	Price
APM956651K	HM525NX Cryostat with UV Disinfection and knife holder	\$22,995

Find out more at thermofisher.com.au