



LabServ Stomacher Bags

The LabServ Stomacher Bags have been specially designed for use with a homogeniser blender. Manufactured to exceed 110µm thickness for additional strength, the LabServ Stomacher Bags provides you with a safe, reliable and economical solution for your homogeniser consumables.

- Gamma sterilised bags
- Robust 110µm construction, very resistant to cuts or punctures
- Reinforced welded seal for extra security and bag integrity
- Safe, reliable and economical solution for sterile homogenising
- Sterility certificate available for complete assurance
- Packaged in individually sealed packs for handling convenience
- Manufactured under ISO9001,2000
- Fits standard 400mL / 3.5L stomacher



Item	Description	Special Price
LBSSTBAG001	Stomacher bags 400mL, sterile, 1000pcs	\$110.00
LBSSTBAG002	Stomacher bags 3500mL, sterile, 400pcs	\$110.00

For more information, contact Thermo Fisher Scientific

© 2018 Thermo Fisher Scientific Inc. Offer valid until 31st December 2018, All rights reserved, Trademarks used are owned as indicated on thermofisher.com.au/trademarks, **1530573633**

Order Placement:

For customer service, call 1300-735-292
 To fax an order, use 1800-067-639
 To email an order, ordersau@thermofisher.com

Contact Us:

For customer service, call 1300-735-292
 For service and calibration, call 1300-736-767
 To order online: thermofisher.com.au